www.swbio.com

Order : order@swbio.com

#11012



NFkB p65(Phospho-Thr435) Antibody

Catalog Number: 11012-1, 11012-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. : Q04206

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl,0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized phosphopeptide derived from human NF κ B p65 around the phosphorylation site of threonine 435 (E-G-TP-L-S).

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatography using non-phosphopeptide corresponding to the phosphorylation site

Specificity/Sensitivity:

NF κ B-p65 (phospho-Thr435) antibody detects endogenous levels of NF κ B-p65 only when phosphorylated at threonine 435

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW: 65kd WB: 1:500~1:1000 IHC: 1:50~1:100 IF:1:100~1:200 Order : order@swbio.com



P-Peptide - +

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue, using NF-кВ p65 (Phospho-Thr435) antibody (#11012).



Peptide

Western blot analysis of extracts from COS7 cells using NF-κB p65 (Ab-435) antibody (#21012, Line 1 and 2) and NF-κB p65 (phospho-Thr435) antibody (#11012, Line 3 and 4)



Immunofluorescence staining of methanol-fixed HeLa cells using NF-κB p65(phospho-Thr435) antibody (#11012,Red)

Background :

NF-kappa-B is a pleiotropic transcription factor present in almost all cell types and is the endpoint of a series of signal transduction events that are initiated by a vast array of stimuli related to many biological processes such as inflammation, immunity, differentiation, cell growth, tumorigenesis and apoptosis. NF-kappa-B is a homo- or heterodimeric complex formed by the Rel-like domain-containing proteins RELA/p65, RELB, NFKB1/p105, NFKB1/p50, REL and NFKB2/p52 and the heterodimeric p65-p50 complex appears to be most abundant one. The dimers bind at kappa-B sites in the DNA of their target genes and the individual dimers have distinct preferences for different kappa-B sites that they can bind with distinguishable affinity and specificity. Different dimer combinations act as transcriptional activators or repressors, respectively. NF-kappa-B is controlled by various mechanisms of post-translational modification and subcellular compartmentalization as well as by interactions with other cofactors or corepressors. NF-kappa-B complexes are held in the cytoplasm in an inactive state complexed with members of the NF-kappa-B inhibitor (I-kappa-B) family. In a conventional activation pathway, I-kappa-B is phosphorylated by I-kappa-B kinases (IKKs) in response to different activators, subsequently degraded thus liberating the active NF-kappa-B complex which translocates to the nucleus. NF-kappa-B heterodimeric p65-p50 and p65-c-Rel complexes are transcriptional activators. The NF-kappa-B p65-p65 complex appears to be involved in invasin-mediated activation of IL-8 expression. The inhibitory effect of I-kappa-B upon NF-kappa-B the cytoplasm is exerted primarily through the interaction with p65. p65 shows a weak DNA-binding site which could contribute directly to DNA binding in the NF-kappa-B complex. Associates with chromatin at the NF-kappa-B promoter region via association with DDX1

References:

Yeh PY, et al. (2004) J Biol Chem. 279(25): 26143-26148. Ryo A, et al. (2003) Mol Cell.12(6): 1413-1426. Baeuerle P A, et al. (1994) Annu Rev Immunol. 12:141-179. Baeuerle P A, et al. (1996) Cell 87:13-20. Haskill S, et al. (1991) Cell 65:1281-1289.