

Catalog Number: 11216-1, 11216-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. :Q14289

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl,0.02% sodium azide and 50% glycerol.

**Storage/Stability:** Store at -20°C/1 year

**Immunogen:** The antiserum was produced against synthesized phosphopeptide derived from Human Pyk2 around the phosphorylation site of tyrosine 402 (D-I-YP-A-E).

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatogramphy using non-phosphopeptide corresponding to the phosphorylation site.

**Specificity/Sensitivity:**Pyk2 (phospho-Tyr402) antibody detects endogenous levels of Pyk2 only when phosphorylated at tyrosine 402.

Reactivity: Human, Mouse, Rat

## Applications:

Predicted MW: 140kd

WB: 1:500~1:1000 IHC: 1:50~1:100 IF:1:100~1:200



РМА — — — +

Peptide - + - -

Western blot analysis of extract from Jurkat cells, untreated or treated with PMA (1ng/ml, 5min), using Pyk2 (Ab-402) antibody (#21209, Lane 1 and 2) and Pyk2 (phospho- Tyr402) antibody (#11216, Lane 3 and 4). www.swbio.com Order : order@swbio.com



Immunohistochemical analysis of paraffin- embedded human breast carcinoma tissue, using Tyk2 (Phospho-Tyr402) antibody (#11216).



Immunofluorescence staining of methanol-fixed HeLa cells Using Tyk2 (Phospho-Tyr402) antibody (#11216).

## Background :

Involved in calcium induced regulation of ion channel and activation of the map kinase signaling pathway. May represent an important signaling intermediate between neuropeptide activated receptors or neurotransmitters that increase calcium flux and the downstream signals that regulate neuronal activity. Interacts with the SH2 domain of Grb2. May phosphorylate the voltage-gated potassium channel protein Kv1.2. Its activation is highly correlated with the stimulation of c-Jun N-terminal kinase activity. Involved in osmotic stress-dependent SNCA 'Tyr-125' phosphorylation

## References:

Gluck SL, et al. (2004) J Clin Invest; 114(12): 1696-1699 Benzing T, et al. (2001) Proc Natl Acad Sci U S A; 98(17): 9784-9789 Tian D, et al. (2002) Mol Cell Biol; 22(8): 2650-2662 Lu Z, et al. (2001) Mol Cell Biol; 21(12): 4016-4031 Krishnan HH, et al. (2006) J Virol; 80(3): 1167-1180