



Tau (Ab-181) Antibody

#21096

Catalog Number: 21096-1, 21096-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. : P10636

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized non-phosphopeptide derived from Human Tau around the phosphorylation site of threonine 181 (P-K-T^P-P-P).


Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity: Tau (Ab- 181) antibody detects endogenous levels of total Tau protein

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW: 48 62 78 kd WB: 1:500~1:1000



Western blot analysis of extract from mouse brain tissue using Tau (Ab-181) antibody (#21096, Lane 1 and 2) and Tau (phospho- Thr181) antibody (#11107, Lane 3 and 4).

Peptide - + - -

P-Peptide - - + -

Background : Promotes microtubule assembly and stability, and might be involved in the establishment and maintenance of neuronal polarity. The C-terminus binds axonal microtubules while the N-terminus binds neural plasma membrane components, suggesting that tau functions as a linker protein between both. Axonal polarity is predetermined by tau localization (in the neuronal cell) in the domain of the cell body defined by the centrosome. The short isoforms allow plasticity of the cytoskeleton whereas the longer isoforms may preferentially play a role in its stabilization.

References: Puig B, et al.(2005) Acta Neuropathol (Berl). 110(3):261-268.

Pérez M, et al.(2005) Neuroscience. 130(2): 339-347.

Fellgiebel A, et al.(2004) Biol Psychiatry. 56(4): 279-283.

del C Alonso A, et al.(2004) J Biol Chem. 279(33): 34873-34881.