www.swbio.com

Order: order@swbio.com





MSK1 (Phospho-Ser376) Antibody



Catalog Number: 11198-1, 11198-2 **Amount:** 50µg/50µl, 100µg/100µl

Swiss-Prot No.: Q75582

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. **Storage/Stability:** Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized phosphopeptide derived from

Human MSK1 around the phosphorylation site of serine 376


Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatogramphy using non-phosphopeptide corresponding to the phosphorylation site.

Specificity/Sensitivity:MSK1 (phospho-Ser376) antibody detects endogenous levels of MSK1 only when phosphorylated at serine 376.

Reactivity: Human, Mouse

Applications:

Predicted MW: 90kd WB: 1:500~1:1000



Western blot analysis of MSK1 phosphorylation expression in HeLa whole cell lysates, The lane on the left is treated with the antigen-specific peptide

Background:

Serine/threonine kinase required for the mitogen or stress-induced phosphorylation of the transcription factors CREB (cAMP response element-binding protein) and ATF1 (activating transcription factor-1). Essential role in the control of RELA transcriptional activity in response to TNF. Directly represses transcription via phosphorylation of 'Ser-1' of histone H2A. Phosphorylates 'Ser-10' of histone H3 in response to mitogenics, stress stimuli and epidemal growth-factor (EGF), which results in the transcriptional activation of several immediate early genes, including proto-oncogenes c-fos/FOS and c-jun/JUN. May also phosphorylate 'Ser-28' of histone H3. Mediates the mitogen- and stress-induced phosphorylation of high mobility group protein 14 (HMG-14).