




HAUSP/USP7

Mouse monoclonal Antibody

#53557

Catalog Number: 53557**Amount:** 100µg/100µl**Swiss-Prot No. :** Q93009**Gene name:** usp7**Gene id:** 7874**Clone Number:** 3E4-D10-E8**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Immunogen:** Purified recombinant human HAUSP / USP7 protein fragments expressed in E.coli.**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of HAUSP / USP7 and does not cross-react with related proteins**Reactivity:** Human**Applications:** Predicted MW: 135kd WB: 1:200

Western blot detection of HAUSP / USP7 in HeLa and K562 cell lysates using HAUSP / USP7 mouse mAb (1:200 diluted), with super ECL. Predicted band size: 135KDa. Observed band size: 135KDa.