



ECFP Tag

Mouse monoclonal Antibody

#T511

Catalog Number: T511

Amount: 100µg/100µl

Swiss-Prot No. :

GENE ID:

Form of Antibody: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Storage/Stability: Store at -20°C/1 year

Purification: affinity-chromatography

Specificity/Sensitivity: This antibody detects endogenous levels of ECFP and ECFP Tag protein


Alternative Names:

Reactivity: ALL

Applications:

Predicted MW:

WB:1:5000-20000



Western blot analysis of recombinant ECFP protein

diluted at 1) 1:5000 2) 1:20000

Background: