


**SMAD3 (Phospho-S425) (1D9)****Mouse monoclonal Antibody****#55231****Catalog Number:** 55231**Amount:** 100µg/100µl**Swiss-Prot No. :**P84022**GENE ID:**4088**Form of Antibody:**Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:**This antibody detects endogenous levels of SMAD3 only when phosphorylated at serine 425.**AlternativeNames:**Mothers against decapentaplegic homolog 3 (MAD homolog 3) (Mad3) (Mothers against DPP homolog 3) (hMAD-3) (JV15-2) (SMAD family member 3) (SMAD 3) (Smad3) (hSMAD3)**Reactivity:** Human,Mouse,Rat**Applications:**

Predicted MW: 52kd

IHC:1:100 -200



Immunohistochemical analysis of paraffin-embedded Human Stomach Carcinoma Tissue using Phospho-Smad3(S425) Mouse mAb

Background: