



Acetyl Histone H3(K9) (2E7)

Mouse monoclonal Antibody


#55121

Catalog Number: 55121**Amount:** 100µg/100µl**Swiss-Prot No. :** P68431**GENE ID:** 8290**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of Acetyl Histone H3(K9) and does not cross-react with related proteins**Alternative Names:** HIST1H3A**Reactivity:** Human, Mouse, Rat

Applications:

Predicted MW: 15kd

IHC: 1:50-200



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl Histone H3 K9 Mouse mAb

Background: