



mStrawberry Tag

Mouse monoclonal Antibody

#T514

**Catalog Number:** T514

**Amount:** 100µg/100µl

**Swiss-Prot No. :**

**GENE ID:**

**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

**Storage/Stability:** Store at -20°C/1 year

**Purification:** affinity-chromatography

**Specificity/Sensitivity:** This antibody detects endogenous levels of mStrawberry and mStrawberry Tag protein

**Alternative Names:**

**Reactivity:** ALL

---

**Applications:**

Predicted MW:                      WB:1:5000



Western blot analysis of recombinant mStrawberry protein

---

**Background:**