



Nano Tag9

Mouse monoclonal Antibody

#T515

Catalog Number: T515

Amount: 100µg/100µl

Swiss-Prot No. :

GENE ID:

Form of Antibody: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Storage/Stability: Store at -20°C/1 year

Purification: affinity-chromatography

Specificity/Sensitivity: This antibody detects endogenous levels of Nano and Nano Tag9 protein

Alternative Names:

Reactivity: ALL

Applications:

Predicted MW:

WB:1:5000



Western blot analysis of recombinant Nano Tag9 protein

Background: