

## Collagen II (7F9)



## Mouse monoclonal Antibody

Catalog Number: 55062 Amount: 100µg/100µl Swiss-Prot No. :P02458

**GENE ID:**1280

Form of Antibody: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM

NaCl) with 0.2% sodium azide, 50%,glycerol **Storage/Stability:** Store at -20°C/1 year **Purification:** affinity-chromatography

Specificity/Sensitivity: This antibody detects endogenous levels of Collagen II and does not corss-react


with related proteins

Alternative Names: Collagen alpha-1(II) chain (Alpha-1 type II collagen) [Cleaved into: Collagen

alpha-1(II) chain; Chondrocalcin] **Reactivity:** Human,Mouse,Rat

Applications:

Predicted MW: 142kd IHC:1:50-200



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma

Tissue using Collagen II Mouse mAb

**Background:**COL2A1 Type II collagen is specific for cartilaginous tissues. It is essential for the normal embryonic development of the skeleton, for linear growth and for the ability of cartilage to resist compressive forces. Belongs to the fibrillar collagen family. Homotrimers of alpha 1(II) chains. 3 isoforms of the human protein are produced by alternative splicing