




Acetyl Lysine

Mouse monoclonal Antibody


#55012

Catalog Number: 55012**Amount:** 100µg/100µl**Swiss-Prot No. :****GENE ID:****Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of Acetyl Lysine and does not cross-react with related proteins**Alternative Names:****Reactivity:** All**Applications:**

Predicted MW: 35kd WB:1:1000-2000 IHC:1:200-500



Western blot analysis of 1) Mouse Brain Tissue, 2) HeLa, 3) HeLa+TSA Treated using Acetyl Lysine Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Acetyl Lysine Monoclonal Antibody

Background: