

-			
Catalog	Number:	24105-1,	24105-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. :Q9UQ07

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl,0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

**Immunogen:** The antiserum was produced against synthesized peptide derived from Human MOK **Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity: MOKantibody detects endogenous levels of total MOK protein


Reactivity: Human

## **Applications:**

Predicted MW: 48kd

WB:1:500~1:2000 IHC

0 IHC:1:100-500



Western blot analysis on HepG2 cell lysate using MOK Antibody

## Background :

MOK a CMGC protein kinase of the RCK family. Possesses a Thr-Glu-Tyr motif in the activation loop domain, as do classical MAP kinases. Associates with and is apparently stabilized by HSP90.