

## APOA1



## Mouse monoclonal Antibody

Catalog Number: 54009 Amount: 100µg/100µl Swiss-Prot No. :P02647

Form of Antibody: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM

NaCl) with 0.2% sodium azide, 50%,glycerol **Storage/Stability:** Store at -20°C/1 year

Immunogen: Purified recombinant human APOA1 protein fragments expressed in E.coli


Purification: affinity-chromatography

Specificity/Sensitivity: This antibody detects endogenous levels of APOA1 and does not corss-react with

related proteins

Reactivity: Human

**Applications:** Predicted MW: 25kd WB: 1:500-2000



Western blot analysis using APOA1 mouse mAb against HEK293T cells transfected with the pCMV6-ENTRY control (1) and pCMV6-ENTRY APOA1 cDNA (2).

## Background:

APOA1: apolipoprotein A-I, it is the major protein component of high density lipoprotein (HDL) in plasma. The protein promotes cholesterol efflux from tissues to the liver for excretion, and it is a cofactor for lecithin cholesterolacyltransferase (LCAT) which is responsible for the formation of most plasma cholesteryl esters. Defects in the Apolipoprotein A I gene are associated with HDL deficiency and Tangier disease. The therapeutic potential of apoA-I has been recently assessed in patients with acute coronary syndromes, using a recombinant form of a naturally occurring variant of apoA-I. The availability of recombinant normal apoA-I should facilitate further investigation into the potential usefulness of apoA-I in preventing atherosclerotic vascular diseases.