

ASS1

Order: order@swbio.com

Mouse monoclonal Antibody

Catalog Number: 54019 Amount: 100µg/100µl Swiss-Prot No.: P00966

Form of Antibody: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM

NaCl) with 0.2% sodium azide, 50%, glycerol Storage/Stability: Store at -20 ℃/1 year

Immunogen: Purified recombinant human ASS1 protein fragments expressed in E.coli


Purification: affinity-chromatography

Specificity/Sensitivity: This antibody detects endogenous levels of ASS1 and does not corss-react with

related proteins

Reactivity: Human, Mouse, Monkey

Applications: Predicted MW: 47kd WB: 1:500-2000 IHC:1:200-1000



Western blot analysis using ASS1 mAb against human ASS1 (AA: 40-236) recombinant protein. (Expected MW is 47 kDa)

Background:

The protein encoded by this gene catalyzes the penultimate step of the arginine biosynthetic pathway. There are approximately 10 to 14 copies of this gene including the pseudogenes scattered across the human genome, among which the one located on chromosome 9 appears to be the only functional gene for argininosuccinate synthetase. Mutations in the chromosome 9 copy of ASS cause citrullinemia. Two transcript variants encoding the same protein have been found for this gene.