



ANKRD1 Antibody

#24137

Catalog Number: 24137-1, 24137-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. :Q15327

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized peptide derived from Human ANKRD1

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity/Sensitivity: ANKRD1 antibody detects endogenous levels of total ANKRD1 protein

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW: 36kd

WB: 1:500~1:2000



Western blot analysis on COLO205 cell lysate using ANKRD1 Antibody

Background :

The protein encoded by this gene is localized to the nucleus of endothelial cells and is induced by IL-1 and TNF-alpha stimulation. Studies in rat cardiomyocytes suggest that this gene functions as a transcription factor.