



#24138

Catalog Number: 24138-1, 24138-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. : Q9H2P0

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized peptide derived from Human ADNP

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity/Sensitivity: ADNP antibody detects endogenous levels of total ADNP protein

Reactivity: Human, Mouse

Applications:

Predicted MW: 124kd

WB: 1:500~1:2000



Western blot analysis on LOVO cell lysate using ADNP Antibody

Background :

ADNP Potential transcription factor. May mediate some of the neuroprotective peptide VIP-associated effects involving normal growth and cancer proliferation.