

ARMX3

Order: order@swbio.com

Antibody

#24139

Catalog Number: 24139-1, 24139-2 Amount: 50μg/50μl, 100μg/100μl Swiss-Prot No. :Q9UH62

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. **Storage/Stability:** Store at -20°C/1 year


Immunogen: The antiserum was produced against synthesized peptide derived from Human ARMX3 **Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity: ARMX3 antibody detects endogenous levels of total ARMX3 protein

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW: 38kd WB:1:500~1:2000 IHC:1:50-200 IF/ICC:1:100-500



Western blot analysis on K562 cell lysate using ARMX3 Antibody

Background:

This gene encodes a member of the ALEX family of proteins which may play a role in tumor suppression. The encoded protein contains a potential N-terminal transmembrane domain and a single Armadillo (arm) repeat. Other proteins containing the arm repeat are involved in development, maintenance of tissue integrity, and tumorigenesis. This gene is closely localized with other family members on the X chromosome. Three transcript variants encoding the same protein have been identified for this gene.