



BBC3 (PUMA) Antibody

#24172

Catalog Number: 24172-1, 24172-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. : Q9BXH1

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized peptide derived from Human BBC3

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity/Sensitivity: BBC3 antibody detects endogenous levels of total BBC3 protein

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW: 20kd

WB: 1:500~1:2000



Western blot analysis of BBC3(PUMA) expression in 293 cells

Background :

Essential mediator of p53/TP53-dependent and p53/TP53-independent apoptosis. Functions by promoting partial unfolding of BCL2L1 and dissociation of BCL2L1 from p53/TP53. Regulates ER stress-induced neuronal apoptosis