



CARD6 Antibody

#24194

Catalog Number: 24194-1, 24194-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. : Q9BX69

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized peptide derived from Human CARD6


Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity: CARD6 Antibody detects endogenous levels of total CARD6

Reactivity: Human

Applications:

Predicted MW: 116kd WB:1:500-2000 IHC:1:50-200 IF/ICC:1:100-500



Western blot analysis on HeLa cell lysate using CARD6 Antibody

Background :

CARD6 May be involved in apoptosis