

CD20



Mouse monoclonal Antibody

Catalog Number: 54034 Amount: 100µg/100µl Swiss-Prot No. :P11836

Form of Antibody: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM

NaCl) with 0.2% sodium azide, 50%,glycerol **Storage/Stability:** Store at -20°C/1 year

Immunogen: Purified recombinant human CD20 protein fragments expressed in E.coli

Purification: affinity-chromatography

Specificity/Sensitivity: This antibody detects endogenous levels of CD20 and does not corss-react with

related proteins

Reactivity: Human

Applications:

Predicted MW: 33kd IHC/ICC:1:200-1000




Figure 1: Immunohistochemical analysis of paraffin-embedded human Tonsil tissues using CD20 mouse mAb.

Background:

CD20 (MS4A1, membrane-spanning 4-domains, subfamily A, member 1) is a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. The CD20 antigen is present on human pre B lymphocytes and on B lymphocytes at all stages of maturation, except on plasma cells. Low level expression of the CD20 antigen has been detected on normal T lymphocytes. The CD20 molecule is involved in regulation of B cell differentiation, presumably via its reported function as a Ca++ channel subunit. And it is known to accelerate the G0 to G1 progression induced by IGF-1.