#54048





Mouse monoclonal Antibody

Catalog Number: 54048

Amount: 100µg/100µl

Swiss-Prot No. : P12830

Form of Antibody:Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

Storage/Stability: Store at -20°C/1 year

Immunogen: Purified recombinant human CDH1 protein fragments expressed in E.coli


Purification: affinity-chromatography

Specificity/Sensitivity: This antibody detects endogenous levels of CDH1 and does not corss-react with related proteins

Reactivity: Human, Mouse, Monkey

Applications:

Predicted MW: 135kd WB:1:500-2000



Western blot analysis using CDH1 mouse mAb against LNCAP (1), A431 (2), DU145 (3), PC-3 (4), PC-12 (5) and T47D(6) cell lysate.

Background: This gene is a classical cadherin from the cadherin superfamily. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein comprised of five extracellular cadherin repeats, a transmembrane region and a highly conserved cytoplasmic tail. Mutations in this gene are correlated with gastric, breast, colorectal, thyroid and ovarian cancer. Loss of function is thought to contribute to progression in cancer by increasing proliferation, invasion, and/or metastasis. The ectodomain of this protein mediates bacterial adhesion to mammalian cells and the cytoplasmic domain is required for internalization. Identified transcript variants arise from mutation at consensus splice sites.