



## FOXP4 Antibody

#24297

**Catalog Number:** 24297-1, 24297-2

**Amount:** 50µg/50µl, 100µg/100µl

**Swiss-Prot No. :** Q8IVH2

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Storage/Stability:** Store at -20°C/1 year

**Immunogen:** The antiserum was produced against synthesized peptide derived from Human FOXP4


**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity/Sensitivity:** FOXP4 Antibody detects endogenous levels of total FOXP4

**Reactivity:** Human, Mouse, Rat

**Applications:**

Predicted MW: 74kd      WB: 1:500-2000



Western blot analysis of extract from JK and HeLa cells using FOXP4

**Background :**

Transcriptional repressor that represses lung-specific expression.