




#24298

Catalog Number: 24298-1, 24298-2**Amount:** 50µg/50µl, 100µg/100µl**Swiss-Prot No. :** Q9H334**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.**Storage/Stability:** Store at -20°C/1 year**Immunogen:** The antiserum was produced against synthesized peptide derived from Human FOXP1**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.**Specificity/Sensitivity:** FOXP1 Antibody detects endogenous levels of total FOXP1**Reactivity:** Human, Mouse, Rat**Applications:**

Predicted MW: 82-90kd WB: 1:500-2000



Western blot analysis of extract from JK, Raji, HeLa and 293 cells using FOXP1 Antibody

Background :

Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium.

Can act with CTBP1 to synergistically repress transcription but CTBP1 is not essential. Essential transcriptional regulator of B cell development