

EGFR (Phospho-Tyr1092) Antibody

#14157

Catalog Number: 14157-1, 14157-2 **Amount:** 50μg/50μl, 100μg/100μl

Swiss-Prot No.: P00533

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. **Storage/Stability:** Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized phosphopeptide derived from Human

EGFR around the phosphorylation site of Tyrosine 1092


Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity:Phospho-EGFR(Tyrosine 1092)Antibody detects endogenous levels of EGFR only when phosphorylated at Tyrosine 1092

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW: 175kd WB:1:500~1:2000 IHC:1:50-200



Western blot analysis of EGFR phosphorylation expression in EGF treated HuvEc whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Background:

EGFR is a receptor tyrosine kinase. Receptor for epidermal growth factor (EGF) and related growth factors including TGF-alpha, amphiregulin, betacellulin, heparin-binding EGF-like growth factor, GP30 and vaccinia virus growth factor. Is involved in the control of cell growth and differentiation.