

## EGFR (Phospho-Tyr1110) Antibody

#14158

**Catalog Number:** 14158-1, 14158-2 **Amount:** 50μg/50μl, 100μg/100μl

Swiss-Prot No.: P00533

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. **Storage/Stability:** Store at -20°C/1 year

**Immunogen:** The antiserum was produced against synthesized phosphopeptide derived from Human

EGFR around the phosphorylation site of Tyrosine 1110

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using


epitope-specific immunogen.

**Specificity/Sensitivity:**Phospho-EGFR( Tyrosine 1110)Antibody detects endogenous levels of EGFR only when phosphorylated at Tyrosine 1110

Reactivity: Human, Mouse, Rat

**Applications:** 

Predicted MW: 175kd WB:1:500~1:2000 IHC:1:50-200



Western blot analysis of EGFR phosphorylation expression in EGF treated A431 whole cell lysates, The lane on the left is treated with the antiqen-specific peptide.

## Background:

EGFR is a receptor tyrosine kinase. Receptor for epidermal growth factor (EGF) and related growth factors including TGF-alpha, amphiregulin, betacellulin, heparin-binding EGF-like growth factor, GP30 and vaccinia virus growth factor. Is involved in the control of cell growth and differentiation.