

## EGFR (Phospho-Ser1071) Antibody

Order: order@swbio.com

#14165

**Catalog Number:** 14165-1, 14165-2 **Amount:** 50μg/50μl, 100μg/100μl

Swiss-Prot No.: P00533

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. **Storage/Stability:** Store at -20°C/1 year

**Immunogen:** The antiserum was produced against synthesized phosphopeptide derived from Human

EGFR around the phosphorylation site of Serine 1071


**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity/Sensitivity:**Phospho-EGFR( Serine 1071)Antibody detects endogenous levels of EGFR only when phosphorylated at Serine 1071

Reactivity: Human, Mouse, Rat

**Applications:** 

Predicted MW: 160kd WB:1:500~1:2000



Western blot analysis of EGFR phosphorylation expression in serum treated HepG2 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

## Background:

EGFR is a receptor tyrosine kinase. Receptor for epidermal growth factor (EGF) and related growth factors including TGF-alpha, amphiregulin, betacellulin, heparin-binding EGF-like growth factor, GP30 and vaccinia virus growth factor. Is involved in the control of cell growth and differentiation.