www.swbio.com

Order: order@swbio.com





PECAM1



Catalog Number: 25601-1, 25601-2 Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No.: P16284

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. Storage/Stability: Store at -20℃/1 year

Immunogen: .Synthesized peptide derived from the Internal region of human PECAM1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using


epitope-specific immunogen.

Specificity/Sensitivity:PECAM1 Antibody detects endogenous levels of total PECAM1

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW:82kd WB:1:500-2000 IHC:1:100-300 IF:1:200-1000



Background: