

Caspase-1 Antibody



Catalog Number: 25530-1, 25530-2 **Amount:** 50µg/50µl, 100µg/100µl

Swiss-Prot No.: P29466

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. **Storage/Stability:** Store at -20°C/1 year

Immunogen: Synthesized peptide derived from the Internal region of human Caspase-1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using


epitope-specific immunogen.

Specificity/Sensitivity: Caspase-1 Antibody detects endogenous levels of total Caspase-1

Reactivity: Human Mouse,Rat

Applications:

Predicted MW:48kd WB:1:500-2000



Western blot analysis of extracts of various cell lines

Background:

This gene encodes a protein which is a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce 2 subunits, large and small, that dimerize to form the active enzyme. This gene was identified by its ability to proteolytically cleave and activate the inactive precursor of interleukin-1, a cytokine involved in the processes such as inflammation, septic shock, and wound healing. This gene has been shown to induce cell apoptosis and may function in various developmental stages. Studies of a similar gene in mouse suggest a role in the pathogenesis of Huntington disease. Alternative splicing results in transcript variants encoding distinct isoforms.