

Bcl-10 Antibody



Catalog Number: 25523-1, 25523-2 Amount: 50μg/50μl, 100μg/100μl Swiss-Prot No. :095999

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. **Storage/Stability:** Store at -20°C/1 year

Immunogen: Synthesized peptide derived from the Internal region of human Bcl-10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using


epitope-specific immunogen.

Specificity/Sensitivity:Bcl-10 Antibody detects endogenous levels of total Bcl-10

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW:36kd WB:1:500-2000 IHC:1:100-300



Western blot analysis of extracts from K562 cells, using BCL10 antibody.

Background:

Involved in adaptive immune response . Promotes apoptosis, pro-caspase-9 maturation and activation of NF-kappa-B via NIK and IKK. May be an adapter protein between upstream TNFR1-TRADD-RIP complex and the downstream NIK-IKK-IKAP complex. Is a substrate for MALT1 .