



EF1 α 1/2 (Acetyl-Lys41) Antibody

#15391

Catalog Number: 15391-1, 15391-2

Amount: 50 μ g/50 μ l, 100 μ g/100 μ l

Swiss-Prot No. : P68104

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized peptide derived from Human EF1 α 1/2 around the acetylation site of Lysine41.


Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity: EF1 α 1/2 (Acetyl-Lys41) Antibody detects endogenous levels of EF1 α 1/2 protein only when acetylated at Lys41.

Reactivity: Human, Mouse, Rat

Applications:

Predicted MW: 50kd WB: 1:500~1:2000



Western Blot analysis of A549 cells using

Acetyl-EF-1 α 1/2 (K41) Antibody

Background :