

P53 (Acetyl-Lys370) Antibody



Catalog Number: 15081-1, 15081-2 **Amount:** 50µg/50µl, 100µg/100µl

Swiss-Prot No.: P04637

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized peptide derived from Human p53 around

the acetylation site of Lysine370.


Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity:p53 (Acetyl-Lys370) Antibody detects endogenous levels of p53 protein only when acetylated at Lys370.

Reactivity: Human **Applications:**

Predicted MW: 53kd

WB:1:500~1:2000 IHC:1:50-100



Western Blot analysis of AD293 cells using

Acetyl-p53 (K370) Antibody

Background :Acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. One of the activated genes is an inhibitor of cyclin-dependent kinases. Apoptosis induction seems to be mediated either by stimulation of BAX and FAS antigen expression, or by repression of Bcl-2 expression. Implicated in Notch signaling cross-over