



**TTR****Mouse monoclonal Antibody****#55259****Catalog Number:** 55259**Amount:** 100µg/100µl**Swiss-Prot No. :**P02766**GENE ID:**7276**Form of Antibody:**Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:**This antibody detects endogenous levels of TTR and does not cross-react with related proteins**AlternativeNames:** Transthyretin (ATTR) (Prealbumin) (TBPA)**Reactivity:** Human**Applications:**

Predicted MW: 16kd WB:1:1000-2000 IHC:1:50-200



Western blot analysis of Human Serum



Immunohistochemical analysis of paraffin-embedded Human Thyroid Tissue

Background: