




Flotillin-1

Mouse monoclonal Antibody

#55258**Catalog Number:** 55258**Amount:** 100µg/100µl**Swiss-Prot No. :** P02751**GENE ID:** 2335**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of Flotillin-1 and does not cross-react with related proteins**AlternativeNames:** FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG**Reactivity:** Human, Mouse, Rat**Applications:**

Predicted MW: 49kd

WB: 1:1000-2000



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain tissue,

Background: