


**MEK2****Mouse monoclonal Antibody****#55266****Catalog Number:** 55266**Amount:** 100µg/100µl**Swiss-Prot No. :** P36507**GENE ID:** 5605**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of MEK2 and does not cross-react with related proteins**Alternative Names:** MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2; MEK 2**Reactivity:** Human, Mouse, Rat**Applications:**

Predicted MW: 45

WB: 1:1000-2000



Western blot analysis of 1) HeLa, 2) 3T3, 3) Rat Brain Tissue

Background: