


**GRP78****Mouse monoclonal Antibody****#55265****Catalog Number:** 55265**Amount:** 100µg/100µl**Swiss-Prot No. :** P11021**GENE ID:** 3309**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of GRP78 and does not cross-react with related proteins**AlternativeNames:** HSPA5; GRP78; 78 kDa glucose-regulated protein; GRP-78; Endoplasmic reticulum luminal Ca(2+)-binding protein grp78; Heat shock 70 kDa protein 5; Immunoglobulin heavy chain-binding protein; BiP**Reactivity:** Human**Applications:**

Predicted MW: 78      WB:1:1000-2000      IHC:1:50-200



Western blot analysis of 1)Hela, 2) Rat LiverTissue

**Background:**